

Taq 2X Master Mix



Catalog #	Size	Concentration	Price	Qty
M0270L	500 reactions (50 µl vol)	2X	\$140.00	<input type="text" value="1"/> 

Categories: [Master Mix Products](#), [Taq DNA Polymerase Products](#)

Applications: [Multiplex PCR](#), [Polymerases for DNA Manipulation](#), [Routine PCR](#) | [More +](#)

Product Information

FAQs & Tech Tips

Protocols & Manuals

Other Tools & Resources

Quality & Safety

Legal Information

- [Description](#)
- [Advantages and Features](#)
- [Properties and Usage](#)
- [Related Products](#)
- [Notes](#)

Description



Already using *Taq*?
Try One *Taq*[®] DNA Polymerase

Taq DNA Polymerase is a thermostable DNA polymerase that possesses a 5'→3' polymerase activity (1,2,3) and

a 5' flap endonuclease activity (4,5).

Taq 2X Master Mix is an optimized ready-to-use solution containing *Taq* DNA Polymerase, dNTPs, MgCl₂, KCl and stabilizers. It is ideally suited to routine PCR applications from templates including pure DNA solutions, bacterial colonies, and cDNA products. It can amplify up to 4 kb from complex genomic DNA or up to 5 kb from lambda DNA.

Highlights

- Robust and reliable reactions
- Tolerates a wide range of templates
- Exceptional value in terms of cost per unit

Product Source

An *E. coli* strain that carries the *Taq* DNA Polymerase gene from *Thermus aquaticus* YT-1

Reagents Supplied

The following reagents are supplied with this product:

	Store at (°C)	Concentration
Magnesium Chloride (MgCl₂) Solution	-20	25 mM

Advantages and Features

Applications

- PCR

- Primer Extension
- Colony PCR

Properties and Usage

Storage Temperature
-20°C

1X Master Mix Composition

10 mM Tris-HCl
50 mM KCl
1.5 mM MgCl₂
0.2 mM dNTPs
5% Glycerol
0.08% IGEPAL® CA-630
0.05% Tween® 20
25 units/ml *Taq* DNA Polymerase
pH 8.6@25°C

Heat Inactivation
No

Unit Assay Conditions

1X ThermoPol® Reaction Buffer, 200 µM dNTPs including [³H]-dTTP and 15 nM primed M13 DNA.

Related Products

Companion Products

- [Deoxynucleotide \(dNTP\) Solution Mix](#)
- [Deoxynucleotide \(dNTP\) Solution Set](#)
- [Standard *Taq* \(Mg-free\) Reaction Buffer Pack](#)
- [Standard *Taq* Reaction Buffer Pack](#)
- [Taq DNA Polymerase with Standard *Taq* Buffer](#)

Materials Sold Separately

- [Magnesium Chloride \(MgCl₂\) Solution](#)

Notes

1. *Taq* 2X Master Mix should be used at a 1X concentration with DNA template and primers in a total reaction volume of 25 or 50 µl.
2. *Taq* 2X Master Mix is stable for fifteen freeze-thaw cycles when stored at -20°C
3. *Taq* 2X Master Mix is also stable for three months at 4°C, so for frequent use, an aliquot may be kept at 4°C.

[FAQs](#)

[Tech Tips](#)

FAQs

1. When should *Taq* DNA Polymerase be used in a primer extension reaction or for PCR?
2. Why is the product a smear when visualized on an agarose gel?
3. Why is there no product when visualized on an agarose gel?
4. The product sequence doesn't completely match the expected sequence. How can this result be improved?
5. Does the presence of Ca²⁺ inhibit PCR?
6. What are the stability and storage requirements of the OneTaq® Master Mixes?

Tech Tips

Did you know most *Taq* reactions amplify more efficiently and robustly when you use a 68°C extension temperature?

[Protocols](#)

[Datacards](#)

Protocols

1. Protocol for *Taq* 2X Master Mix (M0270)

Datacards

The Product Summary Sheet, or Data Card, includes details for how to use the product, as well as details of its formulation and quality controls. The following file naming structure is used to name the majority of these document files: [Catalog Number]Datasheet-Lot[Lot Number]. For those product lots not listed below, please contact NEB at info@neb.com or fill out the [Technical Support Form](#) for appropriate document.

 [M0270Datasheet-Lot0181206](#)

 [M0270Datasheet-Lot0191212](#)

 [M0270Datasheet-Lot0201306](#)

 [M0270Datasheet-Lot0221406](#)

- [Selection Charts](#)
- [Usage Guidelines & Tips](#)
- [Troubleshooting Guides](#)
- [Interactive Tools](#)

Selection Charts

- [DNA Polymerase Selection Chart](#)
- [PCR Selection Chart](#)
- [Thermophilic DNA Polymerases](#)

Usage Guidelines & Tips

- [Activity of Restriction Enzymes in PCR Buffers](#)
- [Guidelines for PCR Optimization with *Taq* DNA Polymerase](#)
- [Guidelines for PCR Optimization with Thermophilic DNA Polymerases](#)

Troubleshooting Guides

- [Taq PCR Kit Troubleshooting Guide](#)
- [PCR Troubleshooting Guide](#)

Interactive Tools

- [NEBioCalculator](#)
- [PCR Selection Tool](#)
- [Tm Calculator](#)

- [Quality Control](#)
- [Safety Data Sheet](#)
- [Datacards](#)

Quality Control

Quality Control Assays

The following Quality Control Tests are performed on each new lot and meet the specifications designated for the product. Individual lot data can be found on the Product Summary Sheet/Datacard or Manual which can be found in the Supporting Documents section of this page. Further information regarding NEB product quality can be found [here](#).

▪ DNase Activity (Labeled Oligo, 3' extension):

The product is tested in a reaction containing a fluorescent labeled double stranded oligonucleotide containing a 3' extension. The percent degradation is determined by capillary electrophoresis.

▪ Endonuclease Activity (Nicking):

The product is tested in a reaction containing a supercoiled DNA substrate. After incubation for 4 hours the percent converted to the nicked form is determined by agarose gel electrophoresis.

▪ PCR Amplification (Master Mix):

The polymerase master mix is tested in a polymerase chain reaction (PCR) using a control template and specific primers,

resulting in the expected product.

Safety Data Sheet

The following is a list of Safety Data Sheet (SDS) that apply to this product to help you use it safely.

 [Taq 2X Master Mix](#)

 [Magnesium Chloride \(MgCl₂\) Solution](#)

Datacards

The Product Summary Sheet, or Data Card, includes details for how to use the product, as well as details of its formulation and quality controls. The following file naming structure is used to name the majority of these document files: [Catalog Number]Datasheet-Lot[Lot Number]. For those product lots not listed below, please contact NEB at info@neb.com or fill out the [Technical Support Form](#) for appropriate document.

 [M0270Datasheet-Lot0181206](#)

 [M0270Datasheet-Lot0191212](#)

 [M0270Datasheet-Lot0201306](#)

 [M0270Datasheet-Lot0221406](#)

 [Legal and Disclaimers](#)

Legal and Disclaimers

This product is covered by one or more patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc (NEB).

While NEB develops and validates its products for various applications, the use of this product may require the buyer to obtain additional third party intellectual property rights for certain applications.

For more information, please contact NEB's Global Business Development team at gbd@neb.com.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

Trademarks

NEW ENGLAND BIOLABS® and THERMOPOL® are registered trademarks of New England Biolabs, Inc.

IGE PAL® is a registered trademark of Rhodia Operations.

TWEEN® is a registered trademark of Uniqema Americas, LLC.